

SYRINGE DOSE PIG

Hoy Exclusive

HOY Scandinavian is proud to present our new product line, which has been developed in order to offer quality products within radiation protection while at the same time being easy on the eyes resulting in a more safe, secure, and pleasant experience by patient as well as staff. All products are produced in anodised aluminium.

SYRINGE DOSE PIG TRANSPORTER

- WITH 6 SYRINGE DOSE PIGS

PART NO: 10.400.01

- Length: 29,5 cm
- Diameter: 10 cm
- Weight: 2,5 kg
- Tungsten shielding 3 mm. pb



ACRYLIC LID FOR TRANSPORTER

PART NO: 10.400.02

- Length: 15.5 cm
- Diameter: 11 cm
- Weight: 385 g



SYRINGE DOSE PIG 1 PIECE

PART NO: 10.400.03

- Length: 15 cm
- Diameter: 3 cm
- Inner diameter: 1.5 cm
- Weight 350 g
- Tungsten shielding 3 mm. pb



LID FOR SYRINGE DOSE PIG

PART NO: 10.400.04

- Length: 9.2 cm
- Diameter: 2.5 cm
- Inner diameter: 3.6 cm
- Weight: 50 g



SYRINGE CARRIER

STAINLESS STEEL WITH 3 MM LEAD

PART NO:10.400.06

- Height: 120 mm
- Diameter: 54 mm
- Weight: 600 g

DISHWASHER PROOF



SHIELDED SYRINGE TRAY

FOR 6 SYRINGES

PART NO:10.500.00

Shielding: 3 mm Pb.
Is stackable.

