

NEW Sliding Syringe Shielding

A new generation of syringe shielding

Hoy Scandinavian is proud to present our new generation of syringe shielding, which we have named "Sliding Syringe Shielding". This has been developed in the same high quality as the other products from the "Hoy Exclusive" Line.

The syringe shield is designed to give the best combination of personal protection and a high functionality, together with fitting any type of needle or infusion valve adjusted to the syringe.



PATENT PENDING!

The sliding syringe shielding is designed with a tungsten shielding of 2,5 mm. and a lead glass window with a density of 5,2 g / cm².

This new syringe shielding are to be separated in two parts, by sliding the sides apart from each other, which hereafter gives the user free access to placing the syringe inside the shielding. It is closed by strong magnets and the syringe can be rotated 360 degrees inside the shielding.

