

## BRITEC Syringe Carriers



Bright Technologies heeft een serie loodkoffertjes voor transport van spuiten ontwikkeld met 3 mm of met 6 mm loodafscherming. De koffertjes bestaan uit een stalen behuizing waarin aan de binnenkant lood is gemonteerd. De binnenzijde bestaat uit een diepgetrokken naadloze kunststof bak. De koffer is verder voorzien van een losse uitneembare binnenbak uit plastic.

De binnenruimte is zo ontworpen dat een 10 ml gevulde spuit met een naald (met kap van 25 mm) hierin geplaatst kan worden.

Het koffertje wordt naar keuze met een vast of met een roterend handvat geleverd. De afstand tussen het deksel en het handvat bedraagt 10 cm.

### Specificaties:

*Behuizing: staal (CR4) met poedercoating*  
*Lood: 3 mm of 6 mm (BSEN12588)*  
*Binnenzijde: plastic ABS*  
*Losse tray: vacuümgevormd plastic HIPS*

Code	Description	Weight	Int. dim.	Ext. dim.
4BT-SCBF3	Small Carrier, 3mm Fixed Handle	4.2 kg	240x55x40 mm <sup>3</sup>	280x95x70 mm <sup>3</sup>
4BT-SCBP3	Small Carrier, 3mm Pivot Handle	4.2 kg		
4BT-SCBF6	Small Carrier, 6mm Fixed Handle	6.0 kg		
4BT-SCBP6	Small Carrier, 6mm Pivot Handle	6.0 kg		
4BT-SCBT	Repl. Tray for Small Carrier Box 5/Box			
4BT-MCBF3	Medium Carrier, 3mm Fixed Handle	4.8 kg	240x75x40 mm <sup>3</sup>	280x120x70 mm <sup>3</sup>
4BT-MCBP3	Medium Carrier, 3mm Pivot Handle	4.8 kg		
4BT-MCBF6	Medium Carrier, 6mm Fixed Handle	7.7 kg		
4BT-MCBP6	Medium Carrier, 6mm Pivot Handle	7.7 kg		
4BT-MCBT	Repl. Tray for Med. Carrier Box 5/Box			
CUSTOM SIZES ARE AVAILABLE ON REQUEST				