



SSC

Shielded Container for Syringe Transport

TECHNICAL SPECIFICATIONS

ENG_Rev.02

The shielded container Mod. SSC is suitable for the transport of shielded syringes.
It is made of lead plate with internal and external stainless steel lining to allow easy decontamination and to warrant a great strength

External dimensions
265 x 95 x 85 (h) mm
h 145 mm with handle up

Internal useful dimensions
200 x 70 x 55 (h) mm

Shielding
6 mm Pb for the two long sides and the bottom
12 mm Pb for the two short sides and the lid

Weight = Kg. 9