



















## RADIATION PROTECTION GLASSES

# RADIATION PROTECTION GLASSES

## STANDARD FRAMES

Standard frames over a range of cost effective designs that provide a comprehensive level of protection.

<p><b>Wraparound</b> Best Fit - Most Faces</p>  <p><b>SUITABLE FOR LARGE PRESCRIPTIONS</b></p> <p>This classic design covers the basics – and then some</p> <p>The lightweight nylon frame features an adjustable bridge that distributes weight evenly over your nose</p> <p>Spring hinges hold the temples snug, and are designed to withstand the rigors of hospital use</p> <p>You'll also love the integrated side shields, which provide solid protection AND excellent peripheral vision</p>	<p><b>Fit Overs</b> Multiple Sizes</p>  <p><b>Available in four sizes</b></p> <p>Fitovers are an oversized wraparound that fits over your existing glasses – or can be worn as a stand-alone splash goggle</p> <p>The top of the frame fits close to your brow for superior splash protection, and side shields are standard</p> <p>Available in S, M, L and XL, Fit Overs accommodate most sizes and styles of glasses</p>	<p><b>Protectors</b> Best Fit - Medium to Large Faces</p>  <p>Whether it's a fluoro procedure in the cath lab or X-rays in the OR, the Protector has you covered</p> <p>The larger lenses provide wrap-around radiation shielding, while the curved frame hugs your face for better splash protection</p> <p>Rubber nose pads and cushioned temple bars provide a secure, comfortable fit during procedures</p>	<p><b>Sport</b> Best Fit - Lenses with Large Rx</p>  <p><b>Available in three sizes</b></p> <p>The Sportcomes with an ultra-soft fixed rubber nose piece for additional comfort</p> <p>A thicker eye frame allows for a wide variety of prescription lenses</p> <p>You can add 0.5mm lead equivalency side shields for additional protection from scatter radiation</p>
<p>COLOUR </p>	<p>COLOUR </p>	<p>COLOUR </p>	<p>COLOUR </p>
<p>ORDER CODE: GLWR</p>	<p>ORDER CODE: GLFT</p>	<p>ORDER CODE: GLPRO</p>	<p>ORDER CODE: GLSW</p>
			<p>WITH SIDE LEAD GLSWSS</p>

<p><b>RG-F10</b> Best Fit - Most Faces</p>  <p>Model F10 utilizes a classic safety style that is solid, sturdy and comfortable</p> <p>The frame comes with perforated brow guards, perforated silicone nose guard and rubberized temple bar for a secure fit</p> <p>The lateral shields protect the orbit of the eye with 0.5mm lead equivalency</p>	<p><b>RG-RK2 Goggles</b> Best Fit - Medium to Large Faces</p>  <p>Model RG-RK2 features vented side ports to reduce fogging, and neoprene foam that protects against splashing and sweat</p> <p>The flame-resistant elastic headband is adjustable to help provide a comfortable fit, whilst the large lens pocket allows for a maximal viewing area</p>	<p><b>RG-Q200</b> Best Fit - Large Faces</p>  <p>Model RG-Q200 offers a stylish wrap-around frame style</p> <p>This frame provides the highest of standards when coming to quality, distortion-free glass, alongside radiation-reducing lenses with 0.75mm lead equivalency</p>	<p><b>375 PX Glasses</b> Best Fit - Medium to Large Faces</p>  <p><b>SUITABLE FOR LARGE PRESCRIPTIONS</b></p> <p>Model 375 offers a lightweight and stylish design</p> <p>The frame is constructed of durable TR-90 nylon for excellent high impact and chemical resistance</p> <p>This model fits best for Small and Medium sized faces</p>
<p>COLOUR </p>	<p>COLOUR </p>	<p>COLOUR </p>	<p>COLOUR </p>
<p>ORDER CODE: GLF10SS</p>	<p>ORDER CODE: GLGOG</p>	<p>ORDER CODE: GLQ200</p>	<p>ORDER CODE: GL375</p>
			<p>WITH SIDE LEAD GL375SS</p>











All our glasses come with plano lenses as standard. Most designs can incorporate a prescription lens. But not all glasses, can take all prescriptions. Please talk to your representative about prescriptions and availability, and please be aware that we can only provide a price once we have received the prescription form. NB, the availability of frames is dependent on the fashion industry and can change quickly and without notice. Please look at the website for an 'up to date' selection.

## RADIATION PROTECTION GLASSES









### LIGHTWEIGHT FRAMES

The lightweight frames are generally designed for smaller faces

Sport Mini		RG-EGM Wraparound		Mini Morpheus		RG-1171 Model	
Best Fit - Smaller Faces		Best Fit - Small and Medium Faces		Best Fit - Smaller Faces		Best Fit - Small and Medium Faces	
							
The Sport comes with an ultra-soft fixed rubber nose piece for additional comfort.  A thicker eye frame allows for a wide variety of prescription lenses.		Model RG-EGM offers a stylish wrap-around frame constructed of high impact TR-90 nylon.  Featuring rubberized soft grip temple bars, saddle bridge nose pad and two-tone temple bars, all for extra comfort.		Bring on the attitude with these rockin' designs!  The Morpheus Mini is the same as our Morpheus but made for smaller faces. Comes with an ultra-soft, fixed rubber nose piece for additional comfort.		Model RG-1171 offers a stylish wrap-around frame with a removable foam insert  The frame is constructed of high impact TR-90 nylon and features rubberized soft grip nose pads and side ventilation ports to reduce fogging.	
COLOUR		COLOUR		COLOUR		COLOUR	
ORDER CODE:	GLMINI	ORDER CODE:	GLEGM	ORDER CODE:	GLMINIMOR	ORDER CODE:	GL1171
WITH SIDE LEAD	GLMINISS			WITH SIDE LEAD	GLMINIMORSS		

### PREMIUM FRAMES

Our premium frames utilise the designs of existing fashion brands. With this comes the increased level of ergonomics and comfort that you would expect from sports frames manufacturers.

Nike Brazen		Nike 7127		Recons		Seals	
Best Fit - Medium to Large Faces		Best Fit - Medium Faces		Best Fit - Smaller to Medium Faces		Best Fit - Small to Medium Faces	
							
The Nike Brazen combines streamlined style with comfort and durability.  Featuring a nylon frame, metal accents and a ventilated nose bridge for reduced fogging.  The secure-wrap temple arms ensure a perfect fit every time.		The 7127 is a lightweight injected frame designed for comfort and durability.  It comes with dual injected temples and channels on the inside of the temples to enhance airflow.		Scout out serious comfort and reliability with the Recon.  The rugged frame is designed to maximize airflow for less fogging and features a soft rubber nose pad for added comfort.		Even under pressure, the Seals hold steady with its anti-slip hydrophilic nose pads and temple tips.  The dependable wrap-around frame comes in stylish colors with bold accents.	
COLOUR		COLOUR		COLOUR		COLOUR	
ORDER CODE:	GLNIKE	ORDER CODE:	GL7127	ORDER CODE:	GLRECO	ORDER CODE:	GLSEAL
WITH SIDE LEAD	GLNIKESS			WITH SIDE LEAD	GLRECOSS	WITH SIDE LEAD	GLSEALSS



Telephone: +44 (0) 1706 830086  
 Fax: +44 (0) 1706 830324  
 Email: sales@rothband.com  
 Web: www.rothband.com

Name: \_\_\_\_\_  
 Order Number: \_\_\_\_\_  
 Hospital: \_\_\_\_\_

## Prescription Glasses Order Form

Please Supply The Following Information

### SINGLE VISION PRESCRIPTION

	SPHERE	CYLINDER	AXIS	PRISM	BASE	PUPILLARY DISTANCE
OD RIGHT EYE						RIGHT
OS LEFT EYE						LEFT

**Information needed:**

1. Distance Rx only (please ensure Rx is current for the year)
2. Pupillary Distance

### BIFOCAL/PROGRESSIVE PRESCRIPTION

	SPHERE	CYLINDER	AXIS	PRISM	BASE	PUPILLARY DISTANCE
OD RIGHT EYE						RIGHT
OS LEFT EYE						LEFT

ADD			<b>Information needed:</b> 1. Distance Rx only (please ensure Rx is current for the year) 2. Pupillary Distance 3. Add Bifocal Power			
-----	--	--	---	--	--	--

### READER ONLY PRESCRIPTION

Add power or magnification: \_\_\_\_\_

**Select type of reader:**

- Bifocal (bottom part of the lens with add power, top part without prescription)
- Single Vision (whole lens will have add power)

Special Instructions (specify if you would like a specific segment height): \_\_\_\_\_