




**NUCLEAR  
SHIELDS**



# COMPACT SOURCE STORAGE SAFES

SC13 - SC14

## Specifications

Finishing material Coated steel  
Shielding material Lead

**Model** **SC1317**  
Shielding thickness 17 mm  
Inside dimensions 334 x 334 x 200 mm (w x d x h)  
Weight 200 Kg

**Model** **SC1432**  
Shielding thickness 32 mm  
Inside dimensions 204 x 204 x 204 mm (w x d x h)  
Weight 196 Kg

**Model** **SC1440**  
Shielding thickness 40 mm  
Inside dimensions 284 x 284 x 360 mm (w x d x h)  
Weight 450 Kg

**Model** **SC1440-2**  
Shielding thickness 40 mm  
Inside dimensions 284 x 284 x 283 mm (w x d x h)  
Weight 420 Kg

**Model** **SC1462**  
Shielding thickness 62 mm  
Inside dimensions 150 x 150 x 150 mm (w x d x h)  
Weight 250 Kg

## Features

- Designed to safely store radioactive materials
- Small and modular concept
- Available in many colors and dimensions
- Lockable with a key
- Can be stacked together or used as a standalone safe
- Bottom frame with legs and/or wheels optional

