

## EDM Premium

EDM Premium is a high-quality ultrasound gel, made in France. EDM Premium is designed to ensure optimal contact between the patient's skin and the transducer.

EDM Premium is a colorless gel with a high viscosity. This gel is pH neutral and has a non-greasy composition. In addition, the EDM Premium is hypoallergenic and non-irritating to the skin.

EDM Premium is supplied in 250 ml semi-rigid bottles. They have an easy to clean snap cap to minimize gel build-up around the opening.

A 5L refill pack is also available.

### Features of the EDM Premium:

- Colorless, transparent, non-sterile Ultrasound Gel
- High Viscosity
- hypoallergenic and non-irritating
- not greasy, does not stain
- pH neutral
- Non-corrosive
- Odorless
- water soluble

