

Britec mini-loodscherf met 1 of 2 loodglasvensters



Deze loodscherfmen hebben een frame uit gepoederd staal, de binnenkant en de bodem zijn afgewerkt met RVS.

Het scherm kan met één of met twee loodglasvensters worden geleverd. Het bovenste loodglas is geplaatst onder een hoek van 45 graden, het optionele tweede venster is gemonteerd in de verticale voorzijde.

Dimensions (mm)	Width	Height	Depth
Base section	275 mm		225 mm
Vertical section	275 mm	265 mm	
Window section	250 mm	200 mm	
Overall height		370 mm	
OTHER SIZES AVAILABLE ON REQUEST			

Code	Description
5BT-BVB6	Mini viewing barrier with 6mm lead shielding and 6mm lead glass window (2mm lead equivalence)
5BT-BVB10	Mini viewing barrier with 10mm lead shielding and 12mm lead glass window (3.5mm lead equivalence)
5BT-BVB10W	Mini viewing barrier with 10mm lead shield and two 12mm lead glass windows, one in cantilevered section and one in vertical front section
5BT-BVBSST	Removable stainless steel tray
5BT-BVBNSS	Mini viewing barrier custom size

