

EDM Krystal

EDM Krystal is a clear, colorless gel with high viscosity and a good value for money. The high viscosity reduces artifacts due to poor contact between the skin and the transducer.

The high viscosity of EDM Krystal Ultrasound Gel ensures even distribution and better contact during procedures. EDM Krystal Ultrasound Gel is distinguished by its homogeneity and colorlessness. This water-soluble product with a non-greasy texture is easy to remove from any surface and is stain-free.

The EDM Krystal bottles have a capacity of 250 ml and have an ergonomic opening at the top that is designed to prevent deposits and simplify dosing during gel application.

Features of the EDM Krystal gel:

- Colorless, transparent, non-sterile Ultrasound Gel
- High Viscosity
- hypoallergenic and non-irritating
- not greasy, does not stain
- pH neutral
- Non-corrosive
- Odorless
- water soluble
- ergonomic cap

