

Medical Imaging



FULL BODY
X-RAY PHANTOM
with real human bones

UNIQUE IN
the world 




Made in  Germany

View product
video here:





X-RAY MANIKIN





Full Body X-Ray Phantom

This model is unique in the world and provides excellent training opportunities for positioning and alignment techniques in projection radiography. It should be part of the basic equipment of any radiographic school. The phantom contains a real human skeleton as well as outlines of larynx, lung, heart and kidneys (organs will create a shadow on the image), which allows taking real X-ray images like in a patient. Using a real skeleton provides even smallest guiding structures which is impossible with a plastic skeleton. During assembly of this phantom we pay special attention to the correct size of joint spaces. All joints are moveably mounted allow positioning in all normal x-ray positions (e.g. frog position, pro- and supination of lower arm). The arms can be moved upwards which makes the phantom suitable for use in all kinds of osseous examinations under CT. Each phantom is hand-made one of a kind; it may differ in size and appearance. Depending on the individual phantom it may have some pathologies, outer shape may differ depending on size of the skeleton. The new version was re-designed in co-operation with a well-known German school for radiographers and fits all needs for education in radiography. This phantom is only sold against proof of medical use. Life size.



Also Partial phantoms available!

Ref.no. 7200

Details:



All Erler-Zimmer phantoms including transport cases



Radiographic positioning doll, plastic skeleton

This model offers all features of model 7200 but includes a plastic skeleton and is due to this only suitable for positioning training. Including a transport and storage case.

Ref.no. 7201
(ohne Abb.)

Details:



Got a broken phantom? Have it repaired at our premises in 4 easy steps:

1. Pick-Up Service

The forwarder comes with an empty transport case and collects your phantom.



2. Repair Cost Estimate

We will provide you with a cost estimate. You decide about a repair then.



3. Repair

The phantom will be repaired fast and reasonable priced at our factory. Result: as **NEW!**



4. Delivery

The "new" phantom will be delivered by the forwarder.

