Headrest, Arm support and Handgrip







General

The head/arm/grip support is especially designed to facilitate a comfortable position for a patient during medical procedures while maintaining proper alignment and preventing collision between the equipment involved and the patient. With the head rest back foam, the head is fixated in supine position during (brain) scans In addition, a hand grip can be used by crossing the hands to optimally fixate the hands and keep the patient confined within the equipment.

Main applications

- Research
- Hospital

Specifications

Overall dimensions 44 x 44 x 33 (LxWxH in cm)



FOR-MED is a development and manufacturing company, specialized in customized medical products for radiology, nuclear medicine and radiotherapy.



