


Galli Mov


- ✓ Training, testing and teaching
- ✓ Preclinical applications
- ✓ Clinical applications

- ✓ Simplify your practice
- Make it easier for patients management
- Encourage the development of new radiopharmaceuticals

Galli Mov is a rolling support designed for IRE ELiT's $^{68}\text{Ge}/^{68}\text{Ga}$ generators. Tightly fitting to the generator base, it allows to move it in the hotcell without having to lift it.





Tray



3 wheels



Item No.	Description	Quantity
1	Tray	1
2	Little wheel	3



Support designed for IRE ELIT's generators only.

SPECIFICATIONS

- **Weight:** 800 g (without generator)
- **Dimensions (W*D*H):** 215 mm * 205 mm * 35 mm
- **Height Galli Mov + generator:** 250 mm
- **Material:** White acetal

PRACTICAL INFORMATION

- **Delivery time:** 4 weeks
- **Information:** radiopharmaceutical@ire-elit.eu