

## ULTRA COVER™

Latex



STERILE

NON  
STERILE

- High durability and elasticity, which makes this cover ideal in the operating room
- Easy adaptation to all types of ultrasound probes

Ref.	Dimensions	Description	Folding	Pack. qty
1-519-3230	3,5 x 20 cm	Cover + 2 elastic bands	Rolled	100
86694*	4 x 30 cm	Cover + 1 adhesive tape	Rolled	10
86675	5 x 30 cm	Cover + 1 adhesive tape	Rolled	10
86667.*	6 x 30 cm	Cover + 1 adhesive tape	Flat	10
86633*	8 x 30 cm	Cover + 1 adhesive tape	Flat	10

Ref.	Dimensions	Description	Folding	Pack. qty
86684	4 x 30 cm	Cover only	Rolled	100
86653.	6 x 30 cm	Cover only	Flat	100

\* Also available with 2 adhesive strips and cable protectors.

## ULTRA COVER™ PIP

Polyisoprene



STERILE

NON  
STERILE

- A polyisoprene cover with elasticity and strength equivalent to a high-quality latex cover

Ref.	Dimensions	Description	Folding	Pack. qty
88694PIP	4 x 30 cm	Cover + cable guard 15 x 150 cm	Rolled	10
88667PIP	6 x 30 cm	Cover + 1 adhesive tape	Flat	10
88633PIP	8 x 30 cm	Cover + 1 adhesive tape	Flat	10
88652PIP	8 x 30 cm	Cover + cable guard 10 x 150 cm	Flat	10

Ref.	Dimensions	Description	Folding	Pack. qty
88663PIP	2.6 x 30 cm	Cover only	Rolled	100
88653PIP	6 x 30 cm	Cover only	Flat	100

## KRYSTAL COVER

Kit Polyethylene



STERILE

- Excellent value for money
- The kit includes: 1 cover, 1 sterile gel packet and 2 elastic bands

Ref.	Dimensions	Description	Material	Folding	Pack. qty
YPE1530WG	15 x 30 cm	Krystal cover kit	PE*	Telescopic	25
YPE1580WG	15 x 80 cm	Krystal cover kit	PE*	Telescopic	25
YPE15125WG	15 x 125 cm	Krystal cover kit	PE*	Telescopic	25
YPE15240WG	15 x 240 cm	Krystal cover kit	PE*	Telescopic	25
YPE1420	14 x 20 cm	Krystal cover	PE*	Telescopic	100
YPE1830	18 x 30 cm	Krystal cover	PE*	Telescopic	100
YPE1880	18 x 80 cm	Krystal cover	PE*	Telescopic	100
YPE18110	18 X 110 cm	Krystal cover	PE*	Telescopic	50
YPE18125	18 x 125 cm	Krystal cover	PE*	Telescopic	50
PEI15125WG	15 x 125 cm	Krystal cover kit with insert card	PE*	Telescopic	24

## KRYSTAL COVER

Kit Polyurethane



- Excellent value for money
- The kit includes: 1 cover, 1 sterile gel packet and 2 elastic bands

Ref.	Dimensions	Description	Folding	Pack. qty
PU1590WG	15 x 90 cm	Krystal cover kit	Telescopic	25
PU1560WG	15 x 60 cm	Krystal cover kit	Telescopic	25
PU1530WG	15 x 30 cm	Krystal cover kit	Telescopic	25
PU15240WG	15 x 240 cm	Krystal cover kit	Telescopic	25
PU15150WG	15 x 150 cm	Krystal cover kit	Telescopic	25
PU15125WG	15 x 125 cm	Krystal cover kit	Telescopic	25
PUI1560WG	15 x 80 cm	Krystal cover kit with insert card	Telescopic	24
PUI15125WG	15 x 125 cm	Krystal cover kit with insert card	Telescopic	24

## POLY COVER™

Kit Polyethylene soft touch



- High-quality ergonomic latex-free cover.
- Very soft to the touch
- Elasticity and strength
- The kit contains: 1 cover, 1 sterile gel packet and 2 elastic bands

Ref.	Dimensions	Description	Folding	Pack. qty
PC1295	10 x 244 cm	Poly Cover™ kit	Telescopic	20
PC1291	13 x 61 cm	Poly Cover™ kit	Telescopic	20
PC1290	15 x 122 cm	Poly Cover™ kit	Telescopic	20
PC1292	15 x 244 cm	Poly Cover™ kit	Telescopic	20

## ADHESION™

Kit Polyethylene soft touch + PRODUCT: adhesive gel



- Cover pre-treated with adhesive gel, eliminating artifacts and residual gel on the probe after use
- Makes clean-up and disinfection fast & easy
- Securely fastens the ultrasound probe cover to the transducer head
- The kit contains: 1 cover, 1 sterile gel packet, 2 elastic bands and 2 adhesive strips

Ref.	Dimensions	Description	Folding	Pack. qty
APC1416	15 x 60 cm	Adhesion™ cover kit	Telescopic	20
APC1290	15 x 120 cm	Adhesion™ cover kit	Telescopic	20
APC1292	15 x 240 cm	Adhesion™ cover kit	Telescopic	20

## PULL UP™

Kit Polyurethane + PRODUCT: insert



- Minimizes cross-contamination with its insert card and telescopic folding
- PullUp™ allows cover to be easily propped open for convenient placement of gel and probe
- PullUp™ handle allows for quick application while maintaining a sterile field
- The kit contains: 1 cover with insert card, 1 sterile gel packet and 2 elastic bands

Ref.	Dimensions	Description	Folding	Pack. qty
1-519-2099OR	14 x 61 cm	Pull Up™ cover kit	Telescopic	24
TBD-2095	14 x 91.5 cm	Pull Up™ cover kit	Telescopic	24
1-519-2097OR	14 x 122 cm	Pull Up™ cover kit	Telescopic	24

## KRYSTAL COVER

Latex



NON STERILE

- The most cost-efficient latex solution for your practice (gynecology, urology, radiology, etc.)

Ref.	Dimensions	Description	Folding	Pack. qty
600249.	3.5 x 20 cm	<b>Non-lubricated cover</b>	Rolled	144
PSE10.	3.5 x 20 cm	<b>Non-lubricated cover</b>	Rolled	144
599249.	3.5 x 20 cm	Lubricated <b>cover</b>	Rolled	144
PSE10R.	3.5 x 20 cm	Lubricated <b>cover</b>	Rolled	144

## ULTRA COVER™

Latex



NON STERILE

- Durable cover adaptable to all types of ultrasound probes

Ref.	Dimensions	Description	Folding	Pack. qty
86697	2.6 x 20 cm	Cover only	Rolled	100
86663	2.6 x 30 cm	Cover only	Rolled	100
650249.	3 x 30 cm	Cover only	Rolled	100

## ULTRA COVER™ PIP

Polyisoprene

Latex-free



NON STERILE

- Polyisoprene, or synthetic latex, eliminates the risk of allergies while maintaining the elasticity and strength of latex

Ref.	Dimensions	Description	Folding	Pack. qty
88663PIP	2.6 x 30 cm	Cover only	Rolled	100
88607PIP	4.4 reduced to 1.8 x 31 cm	Cover only	Rolled	100
670249.	3 x 30 cm	Cover only	Rolled	100

# Echo Clean

*Intermediate-level disinfectant wipes and spray for medical devices*



**In just 5 minutes!**

- Protocol of use adapted to the recommendations for Endocavitary Ultrasound Probes (SF2H 2019).
- Active on Coronavirus, Papillomavirus (HPV-18 and HPV-16)

## KRYSTAL COVER™

Latex



- Ideal for invasive procedures.

Ref.	Dimensions	Description	Folding	Pack. qty
YLA3520	3.5 x 20 cm	Cover only	Rolled	100

## PRO COVER™

Latex



- Ideal for invasive procedures.

Ref.	Dimensions	Description	Folding	Pack. qty
1-519-3230	3.5 x 20 cm	Cover + 2 elastic bands	Rolled	100

## ULTRA COVER™

Latex



- High durability, which makes this cover ideal for the operating room
- Easily adapts to all types of ultrasound probes

Ref.	Dimensions	Description	Folding	Pack. qty
86673.	3 x 30 cm	Cover + 1 adhesive band	Rolled	10
86632.	3.6 x 30 cm	Cover + 1 adhesive band	Rolled	10
86609	4.4 reduced to 1.8 x 31 cm	Cover + 1 adhesive band	Semi rolled	10
86608	4.4 reduced to 1.8 x 31 cm	Cover + 1 adhesive band with a cable cover in polyethylene 8 x 230 cm	Semi rolled	10
86689.	5 reduced to 2 x 36 cm	Cover + 1 adhesive band	Semi rolled	10

## ULTRA COVER™ PIP

Polyisoprene Latex-free



- Elasticity and strength equivalent to high quality latex protection

Ref.	Dimensions	Description	Folding	Pack. qty.
88673PIP	3 x 30 cm	Cover only	Rolled	10
88632PIP	3.6 x 30 cm	Cover only	Rolled	10
88681PIP	3.6 x 30 cm	Cover + 1 adhesive band with a cable cover in polyethylene 8 x 200 cm	Rolled	10

## PULL UP™ SLIMFIT

Kit Polyurethane Latex-free + PRODUCT: insert



- Its cardboard insert guarantees a totally sterile procedure (recommended for IVF)
- The kit contains: 1 protection with insert, 1 sachet of sterile gel and 2 latex free elastics

Ref.	Dimensions	Description	Folding	Pack. qty
1-519-2524	11.9 reduced to 3,8 x 61 cm	Pull Up™ Slimfit protection kit	Telescopic	24
1-519-2548	11.9 reduced to 3,8 x 122 cm	Pull Up™ Slimfit protection kit	Telescopic	24

## Transesophageal probe covers

### ULTRA COVER™ ETO

Kit Latex



- Individually or in a kit, this cover fits all types of ultrasound transesophageal probes
- The kit contains: the cover with a mouthpiece, an Aquasonic gel syringe and a fixation clip.

Ref.	Dimensions	Description	Folding	Pack qty.
86590	2.5 reduced to 1.1 x 100 cm	Cover only	Flat	10
86592	2.5 reduced to 1.1 x 100 cm	Adult kit	Flat	12
86593	1.7 reduced to 1.1 x 65 cm	Pediatric kit	Flat	12

### ULTRA COVER™ PIP ETO

Kit Polyurethane Latex-free



This cover is suitable for all types of transesophageal ultrasound transducers. Available pre-installed on the dispenser (ref. 87111) to avoid setup manipulations, thus limiting the risk of asepsis. The kit contains: the cover with a mouthpiece, an Aquasonic gel syringe and a clip for fixation.



Ref.	Dimensions	Description	Folding	Pack qty.
87106	2.5 reduced to 1.1 x 100 cm	Adult kit	Flat	12
87111	3 reduced to 1.1 x 100 cm	Adult kit - Cover pre-assembled with dispenser	Telescopic	12

### KRYSTAL COVER ETO

Kit Latex



- This kit is suitable for all types of transesophageal ultrasound transducers.
- The cover is pre-assembled with the dispenser to avoid setup manipulations, thus limiting the risk of asepsis.
- The kit contains the cover, the gel, a clip for fixation and a mouthpiece with its elastic band.

Ref.	Description	Folding	Pack qty.
TEEPUKIT	Krystal cover ETO	Flat	12

## Per operative probe cover

### LAPAROSCOPY

Polyethylene Latex-free



- This cover features a tip that adapts to a 10 mm probe.

Ref.	Dimensions	Description	Folding	Pack qty.
16565	5.5 x 400 cm	Cover with polyethylene tip and polyethylene cable guard 10 x 400 cm - requires a 11 mm trocar	Flat + telescopic	5